



ResMed

Actual Size: 140 mm (Width)
x 175 mm (Height)

Ruler is for scaling purposes.


How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirTouch™ F20


Full face mask
Full face mask **for Her**

© 2022 ResMed. All Rights Reserved.
Specifications may change without notice.
Product availability may vary across regions.
1019788/2 2022-07



AirTouch™ F20

Full face mask
Full face mask **for Her**



© 2022 ResMed. All Rights Reserved. Specifications may change without notice.
Product availability may vary across regions.



ResMed

Sample size check

10 mm	20 mm	30 mm	40 mm	50 mm	60 mm	70 mm	80 mm	90 mm	100 mm	110 mm	120 mm	130 mm	140 mm	150 mm	160 mm	170 mm
-------	-------	-------	-------	-------	-------	-------	-------	-------	--------	--------	--------	--------	--------	--------	--------	--------